

## Taste-Related Amino Acid Analysis Using Phenylthiocarbamyl (PTC) in Scallop

Hiroko Seki\*

\*Department of Advanced Food Sciences, College of Agriculture, Tamagawa University, Tokyo, Japan.

Published September 22, 2020

### ABSTRACT

Scallops are consumed globally due to their good and distinctive taste. Particularly, scallops are popular and have high amounts of taste-related amino acids. The amino acid profile should thus be investigated as a significant amount of taste-related amino acids may be found in scallops. Currently, a specialized HPLC with a fluorescent detector is used for amino acid analysis but this method has its limitations due to its low flexibility and high operational cost. In this study, scallop-derived amino acids were derivatized to PTC-amino acids using phenylisothiocyanate (PITC) and analyzed for the taste-related amino acids using a versatile HPLC with a UV detector. I focused on glutamine acid, glycine, alanine, and arginine, as these are the major taste components in scallops. First, I examined 0.5 mM of each, as well as a mixed amino acid solution, and confirmed their detection times. Next, I examined scallop amino acid samples and compared them with standard solutions. For the PTC-amino acid derivatization, 70  $\mu$ L of ethanol, 20  $\mu$ L of trimethylamine, and 20  $\mu$ L of PITC were added to 40  $\mu$ L of sample solution and reacted for 30 min. Then, I added 500  $\mu$ L of 50 mM acetic acid buffer (pH = 6.0):acetonitrile (97:3) to the reaction mixture, filtered using a 0.22- $\mu$ m filter, and analyzed these samples using HPLC. The confirmed detection times for each sample are as follows:  $4.2 \pm 0.2$  min (Glutamine acid),  $7.2 \pm 0.2$  min (glycine),  $8.4 \pm 0.2$  min (alanine), and  $7.8 \pm 0.2$  min (arginine) both in separate and mixture solutions. These detection times were confirmed in both scallop samples. Therefore, I established an approach to determine taste-related amino acids in scallops using HPLC with a UV detector and derivatization to PTC-amino acids.

**Keywords:** Taste-related amino acid, Scallop, HPLC, UV detector, PTC-amino acids

**Corresponding author:** Hiroko Seki, Department of Advanced Food Sciences, College of Agriculture, Tamagawa University, Tokyo, Japan, E-mail: sekhiroko@agr.tamagawa.ac.jp

**Citation:** Seki H. (2020) Taste-Related Amino Acid Analysis Using Phenylthiocarbamyl (PTC) in Scallop. Food Nutr Current Res, 3(S1): 5.

**Copyright:** ©2020 Seki H. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.